


Designed for various manufacturers' electric scalpel, and also for an endoscopic light source.
Brightens you up, by effectively casting light deep into the operating field.




Regulatory reference number: 13B2X00217000001 (light guide)

Light Guide for Electric Scalpels

LG20F

Optical Fiber Model



LG20F-K

For electric scalpels manufactured by CONMED



LG20F-J


For harmonic scalpels

Easy-to-use features to effectively illuminate deep inside an operating field

Light Guide for Electric Scalpels LG20F LG20F Optical Fiber Model

- ◆ Its simple structure allows ease of use and makes its cleaning/disinfection easier.
- ◆ We offer a made-to-order fiber fixation ring for any of electric scalpels produced by several different manufacturers.
- ◆ The model for an endoscopic light source* can be used with the source installed in an operation room.

*Can be attached to the endoscopic light guide adaptor manufactured by OLYMPUS.




Covers a wide range of applications! Light Guide for Electric Scalpels LG20F System


Electric scalpel (CONMED)



LG20F-K



Tip-ring



Electric scalpel

*Inform us of your specific electric scalpel model.

Electric scalpel (MEGADYNE)



LG20F-M




Tip-ring



Harmonic scalpel

Harmonic scalpel



LG20F-J



Tip-ring

Type M

For the light source of optical fiber



Light source for optical fiber




Type O

For an endoscopic light source



Endoscopic light source



Can be connected to the endoscopic light guide adaptor manufactured by OLYMPUS.

■ Models

Light Guide LG20F-K for CONMED Electric Scalpels
LG20F-KM for optical fiber light source
LG20F-KO for endoscopic light source*
Accessories: two rings for the tip of scalpel; Velcro

Light Guide LG20F-M for MEGADYNE Electric Scalpels
LG20F-MM for optical fiber light source
LG20F-MO for endoscopic light source*
Accessories: two rings for the tip of scalpel; Velcro

Light Guide LG20F-J for MEGADYNE Electric Scalpels
LG20F-JM for optical fiber light source
LG20F-JO for endoscopic light source*
Accessories: two rings for the tip of scalpel; Velcro

*Can be connected to the endoscopic light guide adaptor manufactured by OLYMPUS.

■ OPTION



Light source for optical fiber

MegaLight 100
150W Halogen Lamp, weight 1.5 kg
Brightness adjustable, lifetime 1000 hours



STAR MEDICAL, INC.

5F Taiyo-Building 5-14-12 Nishinipponi Arakawa-ku Tokyo 116-0013 Japan TEL:+81-3-5604-1641 FAX:+81-3-5604-1645
E-mail:starmedical@starmedical.co.jp http://www.starmedical.co.jp