

Protective Jacket

For Midget pigs • Rabbits • Cats • Ferrets • Rats

In experiments involving chronic implantation, valuable devices are protected. In long term chronic implantation experiments on both small and large animals, wire leads and connectors drawn out from the body are protected from the animal.

This jacket is made of durable material three dimensionally sewn, and, to prevent shifting of the material, there are reinforcements in several places in the backbone portion. The size is for the upper part of the body, avoids placing undue stress on the animal, and in a natural condition observations can be carried out. The device is placed in the pocket to protect it and further fixed in position with a zip fastener and belt.

For Dogs




Feature

There are a well ventilated mesh type and a sleeve type for use when fastening an electrocardiogram electrode.

Material: Sailcloth (strength)

- With large pocket
Protects a comparatively large device such as an infuser.
Model PJP-D01~10
- With sleeve Height : about 3cm
Handy for use with an ECG electrode.
Model PJS-D01~10
- mesh 1.2 times the price of the standard type.
Model PJM-D01~10

A ... Neck
B ... Chest
C ... Body
D ... Length
E ... Paw



号数	Model	Neck(A)	Chest(B)	Body(C)	Length(D)	Pocket size(E)
1号	PJ-D01	24~27(3)	40~45(5)	40~46(6)	27	5×23
2号	PJ-D02	29~32(3)	46~51(5.5)	45~52(6)	28	6×24
3号	PJ-D03	34~37(3)	52~58(6)	50~57(7)	29	7×25
4号	PJ-D04	39~43(4)	58~65(7)	55~63(8)	31	9×26
5号	PJ-D05	44~49(5)	64~72(8)	60~70(10)	33	10×27
6号	PJ-D06	49~54(5)	70~79(9)	65~76(11)	35	11×30
7号	PJ-D07	54~59(6)	76~86(10)	70~82(12)	37	12×32
8号	PJ-D08	59~65(6)	82~92(10)	75~88(13)	39	13×34
9号	PJ-D09	64~71(7)	88~99(11)	80~94(14)	41	14×36
10号	PJ-D10	69~76(7)	94~106(12)	85~100(15)	43	16×38
Suspending suit		PJ-SS03				

For Rabbits



Feature

Allows the neck to be moved freely. Because it is designed to prevent damage due to "earpicking" movements of the hind legs, that part is particularly strengthened.

号数	Model	Neck	Pocket size
1	PJ-B01		
3	PJ-B03	21~24	7×10
5	PJ-B05		

For Ferrets



Feature

A mesh type, adaptable to gentle movements. There is reinforcement in the back to prevent shifting.

号数	Model	Neck	Pocket size
1	PJ-F01		
3	PJ-F03	19~21	4×3
5	PJ-F05		

For Cats



号数	Model	Neck	Pocket size
1	PC-C01		
3	PJ-C03	23~26	8×12
5	PJ-C05		

For Midget Pigs



Feature

A little tension is applied with the elastic material to prevent the forelegs from coming out.

号数	Model	Neck	Pocket size
1	PJ-MP01		
3	PJ-MP03	54~60	14×20
5	PJ-MP05		

For Rats



Feature

In the backplate beryllium copper plate is used, to which an acrylic resin cylindrical pocket is screwed. An acrylic resin lid is screwed onto it. The thin pipe is lightweight and made of rustproof stainless steel. The pipe has protecting wire fixed by a one-touch operation.



号数	Model	Neck	Pocket size
1	PJ-R01	2	19mmφ×12.5mm
3	PJ-R03	2.5	24mmφ×15mm
5	PJ-R05	3	26mmφ×16.5mm